

CUSTOMER SERVICE BULLETIN




513 Mill Street
Marion, MA 02738-1549
(508) 748-0103 Fax: (508) 748-2029

Date: January 10, 2023
CSB No.: 2023-01
Title: QUIKLOOK FS Pro Software Debut and QUIKLOOK 3.5-FS Hardware


CSB Description:

Teledyne LeCroy Test Services is pleased to announce the latest revision of QUIKLOOK acquisition hardware and software for AOV, MOV, CV, and SOV testing. The new **QUIKLOOK 3.5-FS** (QL3.5-FS) acquisition system, p/n 160800, is a direct replacement for the industry standard QUIKLOOK 3-FS, p/n 160600E. The QL3.5-FS has several advantages over the previous acquisition system; it runs the Windows 10 operating system, has numerous hardware enhancements, and has a completely new version of QUIKLOOK software, **QUIKLOOK FS Pro**.


The user experience is very similar for the QUIKLOOK FS Pro software as compared to the previous version. Menus are the same as the previous version but now appear on a touch-friendly ribbon bar. The replay mode has been enhanced to give the user more flexibility when reviewing and analyzing traces. Gone is the desktop version and desktop Configurator. Now one version exists and the configuration screen for “desktop” looks the same as the “Test” version of the configurator.



QUIKLOOK FS Pro Ribbon Menu



QUIKLOOK FS Pro Configuration Screen



QUIKLOOK FS Pro FlowScanner Tree

CUSTOMER SERVICE BULLETIN



The QUIKLOOK FS Pro version has many new features. The software features zooming enhancements, touch friendly controls, and better graphics. The QUIKLOOK 3.5-FS hardware comes with a built-in keyboard with touchpad, an upgraded touch screen, an additional USB port (3 total), and faster hardware (e.g. internal RAM, motherboard etc). A full list of features for both the QUIKLOOK FS Pro and QUIKLOOK 3.5-FS is available from TLTS personnel upon request. This new software will only run-on networks or on QUIKLOOK 3.5-FS systems (p/n 160800). TLTS recommends updating server installations of QUIKLOOK FS with QUIKLOOK FS Pro since it is backwards compatible with legacy tests as long as those tests contain a .CBD test file. To obtain QUIKLOOK FS Pro a \$0.00 safety related purchase order is required requesting the software.

Customer Action Required:

For technical questions regarding this software or hardware, contact Matthew Ferreira by email (Matthew.Ferreira@Teledyne.com).

Matthew Ferreira

Engineering Manager

Email: Matthew.Ferreira@Teledyne.com

Phone: 508-538-8387

513 Mill Street
Marion, MA 02738-1549
(508) 748-0103 Fax: (508) 748-2029